

Transferring SER-CAT data via Globus

- If you do not have it already, create your Globus end point by installing Globus Connect Personal:
<https://www.globus.org/globus-connect-personal>
- Go to: www.globus.org
- Click on "Login" at the top right section.
- On the next page, select "Argonne National Laboratory" in the organizational login and click 'Continue'. Please login using your Argonne "b#####" or "d#####" number login as you do for your APS user trainings etc. You can contact Argonne IT help desk (help@anl.gov, phone: 630-252-9999) in case of any questions related to your Argonne user login.
- On the next page under "File Manager", choose/type "**APS:DM:SERCAT**" under Collection in the top left box. It may say 'Authentication/Consent required' at the bottom of the page. In that case, click Continue and give the consent for your Argonne identity/login.
- You should be able to see a list of directories with institution names (e.g. UGA, NIH etc.). Click on the corresponding directory with your institution name and your data directories should be listed as sub-directories. Select the data you would like to transfer. In the top right boxes, choose/type your Globus end point and directory path on your local system.
- Start the Globus data transfer by clicking "Start ->".
- **The data should be available for Globus transfer the next day after your collection.**